



University of California
San Francisco

Bone Quality by QCT and HR-pQCT: Translation to Multicenter Clinical Research

Serena Bonaretti, PhD
Department of Radiology and Biomedical Imaging

serenabonaretti.weebly.com


Introduction

- QCT and HR-pQCT are X-ray based imaging techniques used to assess bone quality at a macro-scale (mm) and micro-scale (μm)
- From images, we calculate bone structure, density and strength
- Bone measurements used to evaluate bone quality, investigate age-, sex, and race-related differences, monitor drug therapies, and assess fracture risk
- Currently, QCT and HR-pQCT lack of some degrees of standardization, which prevents them to be used clinically
 - Cross-scanner calibration (QCT)
 - Acquisition procedures (HR-pQCT)

Bone Quality by QCT and HR-pQCT: Translation to Multicenter Clinical Research Scanner Cross-Calibration

Quantitative Computed Tomography (QCT)


- Pad phantom positioned under the patient allows transformation from Hounsfield Units to bone mineral density ($\text{mg}/\text{cm}^3\text{HA}$)



- Volumetric quantification of bone mineral density and geometry separately for cortical and trabecular bone
- Assessment of bone strength through finite element models (geometry and mechanical properties directly derived from images)

Adams 2009; Lang 2010


Why Do We Need Cross-Scanner Calibration?



- Inter-scanner acquisitions
 - Multicenter studies, longitudinal studies with scanner substitution
- Inter-scanner differences are systematic (hardware and software) and patient-specific (beam hardening effect)

Beam Hardening Effect and Patient's Size


- Attenuation of X-rays depends on energy of emitted photons, body tissue and size, and location of the region of interest



Hsieh 2009; Levi 1982

Aim of Cross-Calibration

- To create **correction equations** that allows us to reduce *inter-scanner* errors to *intra-scanner* precision errors
- Different anthropomorphic phantoms are used to cross-calibrate
- Equations are calculated regressing bone mineral density from images acquired on different scanners
- Inter-scanner differences remain and body size has an effect
- There is no current standard solution for QCT scanner cross-calibration




QRM European Spine Phantom

Birnbaum 2014, Carpenter 2014

Anthropomorphic Hip Phantom (AHP)

- Simulates the anatomy and the beam-hardening environment of the pelvis
- Calibration hip to calculate cross-calibration regression
- Test hip inserts to assess the quality of the regression for old and young subjects
- Three body circumferences simulating normal, medium-sized and obese subjects




Experiment

- Anthropomorphic hip phantom scanned on 4 scanners
 - GE VCT 64 System and Siemens Biograph/Definition Flash at UCSF and Mayo Clinic
- Different configuration of age (old/young) and size (no/small/large)
- Cross-calibration procedure
 - Calculated regression line on BMD of phantom *calibration* hip
 - Tested regression line on BMD of phantom *test* hip
 - Applied regression line to BMD of 16 women (age: 64 ± 3 yrs) scanned on GE and Siemens scanners at UCSF


Calculating BMD For Phantom Images

- Same region segmented in all images using affine registration
- One reference image semi-automatically segmented avoiding partial volume effect voxels




Calculating BMD and Body Size For Subjects

- Calculation of BMD (threshold-driven region growing)




Lang 2010, Carpenter 2014

- Calculation of body size (morphologic operators and Hough transform)




Calculating Inter-Scanner Corrections


- Regression lines for BMD of images from different scanners



Calibration Hip - Siemens



Calibration Hip - GE



Regression line
applied to
• test hip
• human subjects


Results – Calibration Hip BMD and Inter-Scanner Corrections

- The most stable scanner with respect to body size was the Siemens at UCSF → reference scanner
- For the other scanners, *intra*-scanner differences of BMD were larger for large body size
- *Inter*-scanner differences of BMD were larger
 - For phantom with large girdle
 - For cortical bone
- In *inter*-scanner regressions, slopes and intercepts were larger for large ring and GE system
- For subjects, we used inter-scanner regression calculated from phantom with old test hip and no girdle

$$\text{BMD}_{\text{GEcrosscal}} = 0.998 \times \text{BMD}_{\text{GE}} - 17.126$$

Results – Testing Corrections on Test Hip and Subjects

- Before correction, inter-scanner differences were significant in all compartments, both for old and young test hip
- After correction, differences reduced considerably, although not significantly
- Mean differences between BMDs decreased after correction, but standard deviation remained the same



	Ct. BMD [mg/cm ³]	Tb. BMD [mg/cm ³]
Before Correction	17.95 ± 17.50	14.42 ± 5.96
After Correction	-5.81 ± 14.35	-3.76 ± 6.67

Discussion

- Novel anthropomorphic hip phantom simulates anatomy and beam hardening effect of the pelvis to calculate cross-scanner corrections
- *Intra-* and *inter*-scanner BMD measurement differences increase when body size increases and are larger for cortical bone
- Cross-scanner correction removed systematic inter-scanner differences, but patient specific differences are still present

Next Steps

- Linear regression does not reduce patient-specific errors in cross-scanner calibration → different approach to the problem
- Currently testing combined cross calibration on 120 patients scanned at UCSF and Mayo Clinic on 4 different scanners

Publications

- **Bonaretti S., Carpenter D.R., Saeed I., Burghardt A.J., Yu L., Bruesewitz M., Khosla S., Lang T. Novel Anthropomorphic Hip Phantom Corrects Systemic Interscanner Differences in Proximal Femoral vBMD.** Phys Med Biol. 59(24), 7819-34. 2014.
- Carpenter R.D., Saeed I., **Bonaretti S., Schreck C., Keyak J.H., Streeper T., Harris T.B., Lang T.F. Inter-scanner Differences in In Vivo QCT Measurements of the Density and Strength of the Proximal Femur Remain After Correction with Anthropomorphic Standardization Phantoms.** Med Eng and Phys. 36(10), 1225-32. 2014.

Acknowledgments


- PI
 - Thomas Lang
- Mayo Clinic
 - L. Yu, M. Bruesewitz, S. Khosla
- Clinical coordinators
 - I. Saeed (UCSF), L. McCready (Mayo Clinic)
- Colleagues
 - D. Carpenter, J. Carballido-Gamio
- Funding
 - NIH/NIAMS 5R01AR060700

Bone Quality by QCT and HR-pQCT: Translation to Multicenter Clinical Research

Standardization of Acquisition – Operator Precision in Reference Line Positioning


High-Resolution Peripheral Quantitative Computed Tomography (HR-pQCT)

- High resolution ($82\mu\text{m}$) *in-vivo* acquisition of distal radius and tibia with low radiation dose ($3\mu\text{Sv}$)
- Morphometric measurements
 - Ct.Th., Ct.Ar., Ct.Po., Tb.N., Tb.Th., Tb.Sp.
- Bone density measurements
 - Tt.BMD, Ct.BMD, Tb.BMD
- Biomechanical measurements
 - L_{failure} , Ct.LF_{prox}, Ct.LF_{dist}




Operator Determines Region to Be Scanned

- Operator positions a reference line at a defined landmark on scout view image
- Scanned region is at a fixed distance from reference line



- Morphometric and densitometric measurements change considerably along bone axis
- Effects on data comparability in cross-sectional studies, in particular multicenter studies
- → **Precise** positioning of reference line



Modified from Boyd, 2008

Previous Literature and Aim


- Scan/rescan total precision errors are smaller
 - For tibia (<5.2%) than for radius (<6.3%)
 - For densitometric (<1.5%) than for structural (<4.5%) measurements
- No study has systematically evaluated the role of HR-pQCT operators in reference line positioning
- To analyze the role of operator ***precision*** in a multicenter HR-pQCT setting
 - Variability of reference line positioning
 - Consequent effects on bone measurements

Boutroy 2005, Khosla 2006, MacNeil 2008, Engelke 2012, Burghardt 2013

Operators and Imaging


- Operators from eight imaging centers

University of Calgary
University of British Columbia
UCSF
Mayo Clinic
Columbia University
Harvard




Université de Lyon
University of Melbourne

- Image acquisition
 - 2 operators at UCSF and Mayo Clinic
 - 56 male and female subjects (71 ± 4 yrs)
 - Double-length stacks (220 slices)
- Simulation outcome measures
 - Reference line positioning
 - Effect on bone parameters



Simulation and Experiment

- Software reproducing acquisition software GUI



- Short-term intra-operator precision
 - Same scout images randomly repeated 3 times
 - 45 images for radius and 48 images for tibia
- Inter-operator precision
 - 50 images for radius and 55 images for tibia

Results - Positioning Precision

POSITIONING PRECISION ERROR	Intra-operator	Inter-operator
	$SD_{RMS}[mm]$	$SD_{RMS}[mm]$
Radius	0.24	0.68
Tibia	0.13	0.30

- Positioning precision error for radius was double than for tibia
- Inter-operator positioning precision was double than inter-operator positioning precision

Results – Bone Parameter Measurements

BONE PARAMETERS ERRORS	Intra-operator $CV_{RMS}[\%]$	Inter-operator $CV_{RMS}[\%]$	Scan-rescan $CV_{RMS}[\%]$
Radius			
Tt.BMD	1.39	3.69	0.89
Ct.Th	3.17	8.40	2.10
Ct.FL _{dist}	2.65	6.37	4.31
Tibia			
Tt.BMD	0.26	0.61	0.43
Ct.Th	0.94	1.97	1.25
Ct.FL _{dist}	0.92	2.02	4.01


- Precision errors for bone parameters followed the trend of positioning precision errors (radius>tibia, inter-op>intra-op)
- Ct.Th had the highest variation for the radius
- Errors due to reference line positioning are equivalent too, or exceed, the in-vivo precision error (scan/rescan)

Discussion and Second Step


- Operator variability in reference line positioning is a major component of in-vivo precision error, for both radius and tibia
- Inter-operator variability is significantly greater than intra-operator variability, for both radius and tibia
- How can we reduce inter-operator variability?

Reference Line – Training and Evaluation

- Standard platform for operator training and certification



- Precise guidelines for reference line positioning



Operators and Experiments

- Six operators involved in the MrOS study with no previous experience with HR-pQCT
- MrOS operators trained on Reference Line – Training and Evaluation and obtained certification
- After certification, operators took part in the same experiments as the experienced operators (intra- and inter-reproducibility)

Results - Positioning Precision

POSITIONING PRECISION ERROR	Intra-operator	Inter-operator
	$SD_{RMS}[mm]$	$SD_{RMS}[mm]$
Radius	0.24/ 0.28	0.68/ 0.30
Tibia	0.13/ 0.11	0.30/ 0.16

Experienced operators/New operators

- Intra-operator errors for new operators were similar errors of experienced operators
- Inter-operator positioning errors for new operators approached intra-operator precision errors

Results – Bone Parameter Measurements

BONE PARAMETERS ERRORS	Intra-operator $CV_{RMS}[\%]$	Inter-operator $CV_{RMS}[\%]$	Scan-rescan $CV_{RMS}[\%]$
Radius			
Tt.BMD	1.39/ 1.50	3.69/ 2.09	0.89
Ct.Th	3.17/ 3.46	8.40/ 4.90	2.10
Ct.FL _{dist}	2.65/ 3.03	6.37/ 3.37	4.31
Tibia			
Tt.BMD	0.26/ 0.31	0.61/ 0.30	0.43
Ct.Th	0.94/ 0.52	1.97/ 1.02	1.25
Ct.FL _{dist}	0.92/ 0.75	2.02/ 1.32	4.01

Experienced operators/[New operators](#)

- Precision errors for bone parameters followed the trend of positioning precision errors (radius>tibia, inter-op~intra-op)

Discussion and Third Step

- Inter-operator positioning precision errors were half for operators who used our training software to learn standard positioning procedures
- Our training platform could be a useful tool to homogenize operator positioning, especially for multicenter studies
- How can we make the tool available for the HR-pQCT community?

Reference Line – Training and Evaluation Web Application

webapps.radiology.ucsf.edu/refline

The screenshot shows a web browser window for the UCSF Reference Line - Training and Evaluation website. The URL in the address bar is webapps.radiology.ucsf.edu/msk/. The page features the UCSF logo and navigation links for About, Documentation, and Contact Us. A large image on the right shows a medical professional in a white coat using a computer system with two monitors, with a patient lying on a scanner table in the foreground. A blue sidebar on the left contains the title "Reference Line Training & Evaluation" and a description: "A tool for the training and evaluation of operator scan positioning for High Resolution Peripheral Quantitative Computed Tomography (HR-pQCT)". At the bottom is a "Launch Web App" button.

Next Steps and Related Works

- Add modules to webapp for hand and knee (OA), and pediatric subjects
- Change radius anatomical landmarks to decrease precision errors
- Anatomically standardized positioning of the region to be scanned (%limb length)
- Detection and correction of patient motion in images, especially radius
- Development of pediatric protocol
- Scanner cross-calibration

Publications

- Bonaretti S., Vilayphiou N., Chan C. M., Yu A., Nishiyama K., Liu D., Boutroy S., Ghasem-Zadeh A., Boyd S.K., Chapurlat R., McKay H., Shane E., Bouxsein M.L., Black D.M., Majumdar S., Orwoll E.S., Lang T.F., Khosla S., Burghardt A.J. **Operator Variability in Scan Positioning Is a Major Component of HR-pQCT Precision Error and Is Reduced by Standardized Training.** Under review in Osteoporosis International.
- Bonaretti S., Holets M., Majumdar S., Lang T.F., Khosla S., Burghardt A.J. **Comparability of HR-pQCT Bone Quality Measures Improved by Scanning Anatomically Standardized Regions.** Ready for submission.

Acknowledgments

- PI
 - Andrew Burghardt
- Experienced HR-pQCT Operators and PIs
 - M. Holets/S. Khosla, K. Nishiyama/E. Shane, D. Liu/H. McKay, S. Boutroy/R. Chapurlat, A. Ghasem-Zadeh/E. Seeman, D. Raymond/S. Boyd, N. Derrico/M. Bouxsein
- MrOS operators
 - D. Cusick, S. Hanson, K. Jacobson, K. Kent, P. Miller, N. Webb
- Webapp
 - Caroline Mai Chan
- Clinical coordinators
 - I. Saeed (UCSF), L. McCready (Mayo Clinic)
- Data management
 - J. Peterson (Mayo Clinic)
- Scanco Medical AG
 - N. Vilayphiou
- Funding
 - NIH/NIAMS 5R01AR060700



University of California
San Francisco

Bone Quality by QCT and HR-pQCT: Translation to Multicenter Clinical Research

Serena Bonaretti, PhD
Department of Radiology and Biomedical Imaging

serenabonaretti.weebly.com